

## Thioaptamer Diagnostic System, Phase II


Completed Technology Project (2010 - 2013)



## Project Introduction

AM Biotechnologies (AM) in partnership with Sandia National Laboratories will develop a Thioaptamer Diagnostic System (TDS) in response to Topic X10.01 Reusable Diagnostic Lab Technology. The TDS will quickly quantify clinically relevant biomarkers in flight using only microliters of virtually any complex sample. The system combines ambient-stable, long-shelf-life affinity agent assays with a handheld microfluidic gel electrophoresis affinity assay quantification technology. The system is easy to use, compatible with operation in microgravity, and designed to permit simultaneous quantification of 32 or more biomarkers from a single astronaut sample. Phase 1 of this project demonstrated that a thioaptamer assay used in the microfluidic instrument can quantify a specific biomarker in serum in the low nanomolar range. AM also successfully identified novel affinity agents to bone specific alkaline phosphatase (BAP) and demonstrated their use to detect BAP using the microfluidic instrument. Phase 2 will expand the number of ambient stable affinity agents and demonstrate a TDS prototype to NASA. AM anticipates that the TDS at the end of Phase 2 will be at TRL 4 to 5. In Phase 3, AM and Sandia will produce flight units for NASA research use on the International Space Station (ISS) as well as for diagnostic use on future long duration missions.

## Primary U.S. Work Locations and Key Partners



Thioaptamer Diagnostic System, Phase II

## Table of Contents

Project Introduction	1
Primary U.S. Work Locations and Key Partners	1
Project Transitions	2
Organizational Responsibility	2
Project Management	2
Technology Maturity (TRL)	2
Technology Areas	3
Target Destinations	3

## Thioaptamer Diagnostic System, Phase II

Completed Technology Project (2010 - 2013)



Organizations Performing Work	Role	Type	Location
AM Biotechnologies, LLC	Lead Organization	Industry	Houston, Texas
● Johnson Space Center(JSC)	Supporting Organization	NASA Center	Houston, Texas

## Primary U.S. Work Locations

Texas

## Project Transitions

**June 2010:** Project Start

## Organizational Responsibility

**Responsible Mission Directorate:**

Space Technology Mission Directorate (STMD)

**Lead Organization:**

AM Biotechnologies, LLC

**Responsible Program:**

Small Business Innovation Research/Small Business Tech Transfer

## Project Management

**Program Director:**

Jason L Kessler

**Program Manager:**


Carlos Torrez

**Principal Investigator:**

Xianbin Yang


## Technology Maturity (TRL)

Start: 4  
 Current: 5  
 Estimated End: 5



## Thioaptamer Diagnostic System, Phase II

Completed Technology Project (2010 - 2013)

**June 2013:** Closed out

**Closeout Summary:** This Phase I project Summary will be formatted as a series of summary statements, with additional comments. The summary statements are intended to focus on the most significant findings and discoveries – things that the project team knows now, but we didn't know at the beginning of the project. 1) Direct Contact Between Gases is an ECLS System enabling capability. The Phase I project focused on CO<sub>2</sub> capture, using direct contact of gases and liquids to facilitate the kinetics of solubility based absorption. At the beginning of the project, we were thinking about the CO<sub>2</sub> capture application, but we were not thinking about ECLS system level impacts. In the course of the phase I effort, we realized that many fundamental unit processes in biological and chemical systems (absorption, evaporation, transpiration, distillation, condensation) need direct contact between a gas phase system and a liquid phase system. Capillary technology that enables gas absorption into a liquid helps all ECLS processes that involve gas absorption – it also enables evaporation, transpiration, distillation, and condensation. Fluid management capability is more broadly applicable than we had originally thought. 2) CO<sub>2</sub> capture using liquid sorbents in microgravity is feasible. The phase I project set out to assess feasibility of the CO<sub>2</sub> capture application. Liquid pumping rates are sufficient, high viscosity is tolerable, contact area is sufficient, flow distribution networks are initially large, but there are design paths to reduce system size. It is feasible to use direct gas/liquid contact for CO<sub>2</sub> capture, and with the right designs and sufficient development, these systems promise to be smaller, more power efficient, and more reliable than current systems. 3) There are other ways to contact gases and liquids – but thin film capillary techniques are new, exciting, and have amazing potential. As a result of this phase I project, we are able to put thin film contactors in context with other fluid management systems. Once placed in context, the power and significance of thin film contact becomes apparent. Two existing methods of gas/fluid management are membrane contactors, and wet spray reactors. Membrane contactors have nearly perfect liquid containment, but unacceptably slow kinetics. Wet spray contactors have great kinetics, but no gas/liquid containment. In the balance between kinetics and containment, thin film contactors are in the "goldilocks zone" – they have sufficient kinetics and sufficient control. This is new and potentially powerful. 4) ECLS system reliability is the key to exploration missions – the key to reliability is having system attributes that favor reliability. Slow moving, gently sweeping, uniform rate, ambient temperature, ambient pressure systems are more likely to operate reliably than systems that involve vacuum systems, high temperatures, hazardous chemicals, and transient operating conditions. If mission planners place a priority on system reliability, ECLS system developers should focus on processes that have the attributes of reliability. 5) The processes with favorable reliability attributes tend to be biological. Biological systems tend to be at ambient temperature, ambient pressure, free of hazardous chemicals, operate at a steady rate, and have the attributes of reliable systems. If Mission Planners want reliable ECLS systems, then ECLS system developers should increase the technical maturity of biological systems. Biological wastewater processes can credibly achieve >98% recovery of water from wastewater. Biological algae based systems can credibly achieve >75% oxygen from CO<sub>2</sub>. 6) The single greatest impact on launch mass of an ECLS system is water. The best way to enable biological water processing is to develop a capillary based method of urine capture that doesn't use pretreat chemicals. The ISS potty uses a rotary phase separator that needs strong pretreat chemicals to prevent biofouling and clogging caused by precipitate build up. These pre-treat chemicals, by design, are intended to kill all microorganisms. Rotary phase separators and biological wastewater systems are

## Technology Areas

## Primary:

- TX06 Human Health, Life Support, and Habitation Systems
  - └ TX06.3 Human Health and Performance
    - └ TX06.3.1 Medical Diagnosis and Prognosis

## Target Destinations

The Sun, Earth, The Moon, Mars, Others Inside the Solar System, Outside the Solar System

## Thioaptamer Diagnostic System, Phase II

Completed Technology Project (2010 - 2013)

